

A-Flex Protective Headgear (Helmet)

USE AND CARE INSTRUCTIONS

The A-Flex Protective Headgear (Helmet) is design to protect patients that have had surgeries to their head.

HOW TO USE YOU'RE A-FLEX PROTECTIVE HEADGEAR

- When applying the A-Flex Protective Headgear, loosen the two side (lateral) velcro straps. Grasp the front and back sections and open or pull out at the same time.
- Place over the patient's head and attach the side velcro straps. Tighten straps. Headgear should be tight enough so it is stable on the patient.
- Under tightening or a loose fit can cause migration of the headgear causing it to tip or come off.

HOW TO CARE FOR YOU'RE A-FLEX PROTECTIVE HEADGEAR

- The A-Flex Protective Headgear (Helmet) can be clean with wiping the inside and outside with Isopropyl Alcohol only. Do not use soap, solvents or perfumes.

IMPORTANT NOTES ABOUT YOUR A-FLEX PROTECTIVE HEADGEAR

Schedule an appointment at our office if any problems with your A-Flex Protective Headgear occur.